

Biossays® C8

Automatic Biochemistry Analyzer

Together with Snibe passion for innovation of **Biossays**[®] C8, the powerful performance and modular scalability enable your lab to realize full automation and maximize efficiency.



www.snibe.com sales@snibe.com





Streamlined workflow

• Decapper module configurable, no manual operation required, largely improve labs efficiency



Super high throughput

- Up to 1600 tests/h (single module)
- Up to 6400 tests/h (four modules combined)



Flexible configuration to meet customized needs

- Capable to link with Laboratory Automation System (TLA/LAS)
- Random combination within 4 units according to testing demand
- Capable to integrate MAGLUMI® X8 (Biolumi® CX8)
- Capable to integrate decapper module



Advanced intelligence



Washing

- 10-step washing technology ensures thorough washing
- Automatic distribution of warm water, acid and alkali liquid reduces the liquid consumption (Assay dependent)
- Independent agitation mechanism of washing liquid to reduce carryover and ensure a more thorough cleaning

Mixing and Incubation

- Full mixing with adjustable speed, the mixing speed of different reagents can be set independently
- Permanent quartz cuvettes, maintenance-free and no replacement required
- Solid heating with stable temperature control at 37±0.2°C



Anemia

Transferrin (TRF) Iron UIBC *G6PD *Ferritin

Pancreatic

α-AMY LIP *P-AMY

Diabetes

GLU
LAC
HbA1c
GSP
GA
β-Hydroxybutyrate
*Adiponectin (ADPN)

Rheumatism

ASO RF

Hepatic
ТВА
ALT (SGPT)
AST (SGOT)
ALP
GGT
TBIL (Vanadate)
DBIL (Vanadate)
TP
ALB
Ammonia (AMM)
PA
CHE
AFU
5´-NT
Haptoglobin (HPT)
TBIL (DPD)
ADA
LAP
GR
*Monoamine Oxidase
*DBIL (DCA)
*CG (Latex)
*CG (CEDIA)
*GLDH
*CER

IgG C3 C4 *C1q *Kappa Light Chain *Lambda Light Chain	a-HBI LDH Hcy LDH1 *D-Di *ACE *IMA *MYC
Renal	*mAS
Cr (CREA) Uric Acid (UA) Urea	Inorg
Cys-C α1-Microglobulin β2-Microglobulin Urine / CSF Protein mALB RBP	Ca (C P (PH Mg Ca (A CO ₂
NAG	
*α2-Macroglobulin	Inflar

Immune

*Urine-IgG

ΙgΜ

IgA

Micro-sampling technology to ensure the sampling accuracy of the low-volume sample (1.5 μL) Liquid level, clot and collision detection technologies reduce sampling error rate

• Reagent probe: liquid level detection, collision detection

Measuring

Pipetting

- The precise control of temperature and voltage of light source lamp guarantee the small fluctuation of the light source and long life circle of lamp
- Concave holographic grating post-splitting technology: high resolution and low stray light
- High-performance measurement channels: high signal-to-noise ratio, small ripple and small fluctuation of temperature

Cardiac	Lipids
CK CK-MB a-HBDH LDH Hcy LDH1 *D-Dimer *ACE *IMA *MYO	HDL-C LDL-C TC TG ApoE ApoA1 ApoB Lp (a) NEFA
*mAST	ISE
*hs-CRP	Na⁺ K⁺
Inorganic Ion	CI
Ca (OCPC) P (PHOS) Mg Ca (Arsenazo III) CO ₂	Ca²⁺ pH
Inflammation	
CRP (Full Range)	

*Available soon

Technical Specifications

General information				
Test throughput (Theoretical max)	Single module: 1600 tests/hour Four modules combined: 6400 tests/hour ISE module: 300 tests/hour			
Sample positions	300			
Sample barcode types	Code128, Code39, Code93, Codabar, 2/5 Interleaved			
Emergency mode	STAT function available Continuous operation: additional sample loading during operation			
QC types	Batch QC, Month QC			
Dimensions	Biochemical module: 133*118*135 cm, 560 kg Sample area: 44*118*99 cm, 90 kg Decapper module: 58*118*135 cm, 170 kg			
Biochemistry module				
Sample types	Serum, Plasma, Urine, CSF			
Sample volume	1.5 - 25 μL in 0.1 μL steps			
Sampling probe	Automatic washing, liquid level detection, clot detection, collision detection			
Reagent compartments	72 positions each for R1 and R2			
Reagent temperature	Working temperature: 8°C - 12°C Storage temperature: 2°C - 8°C			
Reagent probe	Automatic washing, liquid detection, collision detection			
Reaction cuvettes	362 cuvettes in total: 181 cuvettes each for two segments			
Analysis time	10 minutes, 22 minutes			
Reaction temperature	$37.0^{\circ}C \pm 0.2^{\circ}C$			
Reaction volume	80 μL - 250 μL			
Wavelength	Photometer: 13 fixed wavelengths (340, 380, 405, 450, 480, 505, 546, 570, 600, 660, 700, 750, 800 nm)			
Light source	12 V, 100 W halogen lamp			
ISE module (configuration available)				
Sample volume	90 μL			
Information recognition	RFID			
Electrode types	K ⁺ , Na ⁺ , Cl ⁻ , iCa ²⁺ , pH, reference electrode			
Decapper module (configuration available)				
Throughput	Decapper speed: 500 tubes/hour Capper speed: 400 tubes/hour Combined speed for capper and decapper: 400 tubes/hour			

Biossays[®], **Biolumi**[®] and **MAGLUMI**[®] are trademarks of Snibe. All other product names and trademarks are the property of their respective owners.

Snibe Diagnostic

Shenzhen New Industries Biomedical Engineering Co., Ltd. (SNIBE Co., Ltd.)

No.23, Jinxiu East Road, Pingshan District, 518122 Shenzhen, P.R. China

Tel: +86 755 26501514 Email: sales@snibe.com



Fax: +86 755 26654850 Web: www.snibe.com